

## BLK 抗原（重组蛋白）

中文名称：BLK 抗原（重组蛋白）

英文名称：BLK Antigen (Recombinant Protein)

别名：MODY11

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 300 amino acids of human BLK

技术规格：

<b>Full name:</b>	B lymphoid tyrosine kinase
<b>Synonyms:</b>	MODY11
<b>Swissprot:</b>	P51451
<b>Gene Accession:</b>	BC007371
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a nonreceptor tyrosine-kinase of the src family of proto-oncogenes that are typically involved in cell proliferation and differentiation. The protein has a role in B-cell receptor signaling and B-cell development. The protein also stimulates insulin synthesis and secretion in response to glucose and enhances the expression of several pancreatic beta-cell transcription factors.