

## CATSPER3 抗原(重组蛋白)

中文名称: CATSPER3 抗原(重组蛋白)

英文名称: CATSPER3 Antigen (Recombinant Protein)

别 名: CACRC

储 存: 冷冻(-20℃)

相关类别: 抗原

概 述:

Fusion protein corresponding to a region derived from 299-398 amino acids of human CATSPER3

## 技术规格:

Full name:	cation channel, sperm associated 3				
Synonyms:	CACRC				
Swissprot:	Q86XQ3				
Gene Accession:	BC110384				
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE				
Expression system:	Escherichia coli				
Tags:	His tag C-Terminus, GST tag N-Terminus				
Background:	CatSper3, is a protein which in humans is encoded by the CATSPER 3 gene. CatSper3 is a member of thecation channels of sperm famil y of proteins. The four proteins in this family together form a Ca2+ -permeant ion channel specific essential for the correct function of s perm cells. Voltage-gated calcium channel that plays a central role in calcium-dependent physiological responses essential for successful				



fertilization,	such	as sperm	hyperactivation,	acrosome	reaction	and	С
hemotaxis t	owards	the ooc	yte.				