

CDK5RAP3 抗原（重组蛋白）

中文名称： CDK5RAP3 抗原（重组蛋白）

英文名称： CDK5RAP3 Antigen (Recombinant Protein)

别名： C53; IC53; LZAP; HSF-27; MST016; PP1553; OK/SW-cl.114

相关类别： 抗原

储存： 冷冻（-20℃）

概述：

Fusion protein corresponding to a region derived from 248-506 amino acids of human CDK5RAP3

技术规格：

Full name:	CDK5 regulatory subunit associated protein 3
Synonyms:	C53; IC53; LZAP; HSF-27; MST016; PP1553; OK/SW-cl.114
Swissprot:	Q96JB5
Gene Accession:	BC009957
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein that has been reported to function in signaling pathways governing transcriptional regulation and cell cycle progression. It may play a role in tumorigenesis and metastasis. A pseudogene of this gene is located on the long arm of chromosome 20. Alternative splicing results in multiple transcript variants that encode different isoforms.

mlbio 梅联生物
Good elisakit producers