

## CENPC 抗原（重组蛋白）

中文名称： CENPC 抗原（重组蛋白）

英文名称： CENPC Antigen (Recombinant Protein)

别名： MIF2; hcp-4; CENP-C; CENPC1

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to C terminal 250 amino acids of human CENPC	
---	--

技术规格：

<b>Full name:</b>	centromere protein C
<b>Synonyms:</b>	MIF2; hcp-4; CENP-C; CENPC1
<b>Swissprot:</b>	Q03188
<b>Gene Accession:</b>	BC030695
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Centromere protein C 1 is a centromere autoantigen and a component of the inner kinetochore plate. The protein is required for maintaining proper kinetochore size and a timely transition to ana phase. A putative pseudogene exists on chromosome 12. Component of the CENPA-NAC (nucleosome-associated) complex, a complex that plays a central role in assembly of kinetochore proteins, mitotic progression and chromosome segregation. The CENPA-NA

C complex recruits the CENPA-CAD (nucleosome distal) complex and may be involved in incorporation of newly synthesized CENPA into centromeres.