

CHDH 抗原（重组蛋白）

中文名称： CHDH 抗原（重组蛋白）

英文名称： CHDH Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概 述：

Fusion protein corresponding to C terminal 250 amino acids of human CHDH

技术规格：

Full name:	choline dehydrogenase
Swissprot:	Q8NE62
Gene Accession:	BC034502
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a choline dehydrogenase that localizes to the mitochondrion. Variations in this gene can affect susceptibility to choline deficiency. A few transcript variants have been found for this gene, but the full-length nature of only one has been characterized to date.