

## CYP2W1 抗原（重组蛋白）

中文名称： CYP2W1 抗原（重组蛋白）

英文名称： CYP2W1 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概 述：

Fusion protein corresponding to a region derived from C terminal 250 amino acids of human CYP2W1

技术规格：

<b>Full name:</b>	cytochrome P450, family 2, subfamily W, polypeptide 1
<b>Swissprot:</b>	Q8TAV3
<b>Gene Accession:</b>	BC025761
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids.