

## DIRAS3 抗原（重组蛋白）

中文名称：DIRAS3 抗原（重组蛋白）

英文名称：DIRAS3 Antigen (Recombinant Protein)

别名：ARHI; NOEY2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-226 amino acids of human DIRAS3

技术规格：

<b>Full name:</b>	DIRAS family, GTP-binding RAS-like 3
<b>Synonyms:</b>	ARHI; NOEY2
<b>Swissprot:</b>	O95661
<b>Gene Accession:</b>	BC005362
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene is a member of the ras superfamily, and is expressed in normal ovarian and breast epithelial cells, but not in ovarian and breast cancers. It is an imprinted gene, with monoallelic expression of the paternal allele, which is associated with growth suppression. Thus, this gene appears to be a putative tumor suppressor gene.

ene whose function is abrogated in ovarian and breast cancers.