

EPHB6 抗原（重组蛋白）

中文名称：EPHB6 抗原（重组蛋白）

英文名称：EPHB6 Antigen (Recombinant Protein)

别名：HEP

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 250 amino acids of human EPHB6

技术规格：

Full name:	EPH receptor B6
Synonyms:	HEP
Swissprot:	O15197
Gene Accession:	BC110606
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of a family of transmembrane proteins that function as receptors for ephrin-B family proteins. Unlike other members of this family, the encoded protein does not contain a functional kinase domain. Activity of this protein can influence cell adhesion and migration. Expression of this gene is downregulated during tumor progression, suggesting that the protein may suppress tumor invasion and metastasis. Alt

ernative splicing results in multiple transcript variants.