

## EZR 抗原（重组蛋白）

中文名称： EZR 抗原（重组蛋白）

英文名称： EZR Antigen (Recombinant Protein)

别名： CVL; CVIL; VIL2; HEL-S-105

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 2-295 amino acids of human EZR

技术规格：

<b>Full name:</b>	ezrin
<b>Synonyms:</b>	CVL; CVIL; VIL2; HEL-S-105
<b>Swissprot:</b>	P15311
<b>Gene Accession:</b>	BC013903
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The cytoplasmic peripheral membrane protein encoded by this gene functions as a protein-tyrosine kinase substrate in microvilli. As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this

gene. Alternatively spliced variants have also been described for this gene.