

HOMER1 抗原（重组蛋白）

中文名称： HOMER1 抗原（重组蛋白）

英文名称： HOMER1 Antigen (Recombinant Protein)

别 名： HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C

储 存： 冷冻（-20℃）

相关类别： 抗原

概 述：

Fusion protein corresponding to C terminal 125 amino acids of human HOMER1

技术规格：

Full name:	homer homolog 1
Synonyms:	HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C
Swissprot:	Q86YM7
Gene Accession:	BC015502
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the homer family of dendritic proteins. Members of this family regulate group 1 metabotropic glutamate receptor function.