

IL1RL1 抗原（重组蛋白）

中文名称：IL1RL1 抗原（重组蛋白）

英文名称：IL1RL1 Antigen (Recombinant Protein)

别名：Interleukin 1 receptor-like 1; T1; ST2; DER4; ST2L; ST2V; FIT-1; IL33R

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 250 amino acids of human IL1RL1

技术规格：

Full name:	Interleukin 1 receptor-like 1
Synonyms:	T1; ST2; DER4; ST2L; ST2V; FIT-1; IL33R
Swissprot:	Q01638
Gene Accession:	BC030975
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the interleukin 1 receptor family. Studies of the similar gene in mouse suggested that this receptor can be induced by proinflammatory stimuli, and may be involved in the function of helper T cells. This gene, interleukin 1 receptor, type I (IL1R1), interleukin 1 receptor, type II (IL1R2) and interleukin 1 receptor-like 2 (IL1RL2) form a cytokine receptor gene cluster in a region mapped to chromosome 2q12. Alternative splicing of this gene results in multiple transcript variants.

mlbio 梅联生物
Good elisakit producers