

IL18RAP 抗原（重组蛋白）

中文名称：IL18RAP 抗原（重组蛋白）

英文名称：IL18RAP Antigen (Recombinant Protein)

别名：Interleukin 18 receptor accessory protein; ACPL; CD218b; IL18RB; CDw218b

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 378-599 amino acids of human IL18RAP

技术规格：

Full name:	Interleukin 18 receptor accessory protein
Synonyms:	ACPL; CD218b; IL18RB; CDw218b
Swissprot:	O95256
Gene Accession:	BC106765
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Interleukin 18 receptor accessory protein, also known as IL18RAP and CDw218b (cluster of differentiationw218b), is a human gene. The protein encoded by this gene is an accessory subunit of the heterodimeric receptor for IL18. This protein enhances the IL18 binding activity of IL18R1 (IL1RRP), a ligand binding subunit of IL18 receptor. The coexpression of IL18R1 and this protein is required for the activation of NF-κB and MAPK8 (JNK) in response to IL18.

mlbio 梅联生物
Good elisakit producers