

MAP4K2 抗原（重组蛋白）

中文名称：MAP4K2 抗原（重组蛋白）

英文名称：MAP4K2 Antigen (Recombinant Protein)

别名：GCK; BL44; RAB8IP

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 250 amino acids of human MAP4K2

技术规格：

Full name:	mitogen-activated protein kinase kinase kinase kinase 2
Synonyms:	GCK; BL44; RAB8IP
Swissprot:	Q12851
Gene Accession:	BC047865
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the serine/threonine protein kinase family. Although this kinase is found in many tissues, its expression in lymphoid follicles is restricted to the cells of germinal centre, where it may participate in B-cell differentiation. This kinase can be activated by TNF-alpha, and has been shown to specifically activate MAP kinases. This kinase is also found to interact with TNF receptor-associated factor 2 (TRAF2), whi

ch is involved in the activation of MAP3K1/MEKK1.