

MLLT3 抗原（重组蛋白）

中文名称：MLLT3 抗原（重组蛋白）

英文名称：MLLT3 Antigen (Recombinant Protein)

别名：MLLT3, super elongation complex subunit; AF9; YEATS3

相关类别：抗原

储存：冷冻（-20℃）

概述：

Fusion protein corresponding to a region derived from 322-568 amino acids of human MLLT3

技术规格：

| | |
|---------------------------|---|
| Full name: | MLLT3, super elongation complex subunit |
| Synonyms: | AF9; YEATS3 |
| Swissprot: | P42568 |
| Gene Accession: | BC036089 |
| Purity: | >85%, as determined by Coomassie blue stained SDS-PAGE |
| Expression system: | Escherichia coli |
| Tags: | His tag C-Terminus, GST tag N-Terminus |
| Background: | Component of the super elongation complex (SEC), a complex required to increase the catalytic rate of RNA polymerase II transcription by suppressing transient pausing by the polymerase at multiple sites along the DNA. |