

MTM1 抗原（重组蛋白）

中文名称：MTM1 抗原（重组蛋白）

英文名称：MTM1 Antigen (Recombinant Protein)

别名：CNM; MTMX; XLMTM

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 240-484 amino acids of human MTM1

技术规格：

Full name:	myotubularin 1
Synonyms:	CNM; MTMX; XLMTM
Swissprot:	Q13496
Gene Accession:	BC030779
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a dual-specificity phosphatase that acts on both phosphotyrosine and phosphoserine. It is required for muscle cell differentiation and mutations in this gene have been identified as being responsible for X-linked myotubular myopathy