

## 兔抗 CRABP2 多克隆抗体

中文名称：兔抗 CRABP2 多克隆抗体

英文名称： Anti-CRABP2 rabbit polyclonal antibody

别名： RBP6; CRABP-II

相关类别： 一抗

抗原： CRABP2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the retinoic acid (RA, a form of vitamin A) binding protein family and lipocalin/cytosolic fatty-acid binding protein family. The protein is a cytosol-to-nuclear shuttling protein, which facilitates RA binding to its cognate receptor complex and transfer to the nucleus. It is involved in the retinoid signaling pathway, and is associated with increased circulating low-density lipoprotein cholesterol. Alternatively spliced transcript variants encoding the same protein have been found for this gene.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	CRABP2
<b>Immunogen:</b>	Fusion protein of human CRABP2
<b>Full name:</b>	cellular retinoic acid binding protein 2
<b>Synonyms :</b>	RBP6; CRABP-II
<b>SwissProt:</b>	P29373
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human breast cancer and human cervical cancer
<b>IHC Recommend dilution:</b>	50-200

