

兔抗 CRB3 多克隆抗体

中文名称：兔抗 CRB3 多克隆抗体

英文名称： Anti-CRB3 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： CRB3

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the Crumbs family of proteins. This gene is widely expressed in epithelial tissues where the encoded protein isoforms play various roles such as the control of cytokinesis and ciliogenesis or the formation of tight junctions. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	CRB3
Immunogen:	Fusion protein of human CRB3
Full name:	crumbs 3, cell polarity complex component
SwissProt:	Q9BUF7
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human tonsil;
IHC Recommend dilution:	25-100

