

兔抗 CREB1 多克隆抗体

中文名称：兔抗 CREB1 多克隆抗体

英文名称：Anti-CREB1 rabbit polyclonal antibody

别名：CREB

相关类别：一抗

抗原：CREB1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a transcription factor that is a member of the leucine zipper family of DNA binding proteins. This protein binds as a homodimer to the cAMP-responsive element, an octameric palindrome. The protein is phosphorylated by several protein kinases, and induces transcription of genes in response to hormonal stimulation of the cAMP pathway. Alternate splicing of this gene results in two transcript variants encoding different isoforms.
Applications:	ELISA, WB
Name of antibody:	CREB1
Immunogen:	Fusion protein of human CREB1

Full name:	cAMP responsive element binding protein 1
Synonyms :	CREB
SwissProt:	P16220
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	37 kDa
WB Positive control:	NIH/3T3 cells
WB Recommended dilution:	200-1000

