

## 兔抗 CREG1 多克隆抗体

中文名称：兔抗 CREG1 多克隆抗体

英文名称：Anti-CREG1 rabbit polyclonal antibody

别名：cellular repressor of E1A stimulated genes 1 CREG

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CREG1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The adenovirus E1A protein both activates and represses gene expression to promote cellular proliferation and inhibit differentiation. The protein encoded by this gene antagonizes transcriptional activation and cellular transformation by E1A. This protein shares limited sequence similarity with E1A and binds both the general transcription factor TBP and the tumor suppressor pRb in vitro. This gene may contribute to the transcriptional control of cell growth and differentiation.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	CREG1

<b>Immunogen:</b>	Synthetic peptide of human CREG1
<b>Full name:</b>	cellular repressor of E1A stimulated genes 1
<b>Synonyms:</b>	CREG
<b>SwissProt:</b>	O75629
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	24 kDa
<b>WB Positive control:</b>	RAW264.7 cell lysate
<b>WB Recommended dilution:</b>	500-2000

