

## 兔抗 CRELD1 多克隆抗体

中文名称：兔抗 CRELD1 多克隆抗体

英文名称：Anti-CRELD1 rabbit polyclonal antibody

别名：AVSD2; CIRIN

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CRELD1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a subfamily of epidermal growth factor-related proteins. The encoded protein is characterized by a cysteine-rich with epidermal growth factor-like domain. This protein may function as a cell adhesion molecule. Mutations in this gene are the cause of atrioventricular septal defect. Alternate splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CRELD1
<b>Immunogen:</b>	Synthetic peptide of human CRELD1

<b>Full name:</b>	cysteine-rich with EGF-like domains 1
<b>Synonyms :</b>	AVSD2; CIRRIN
<b>SwissProt:</b>	Q96HD1
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human brain
<b>IHC Recommend dilution:</b>	50-200

