

ORM1 抗原（重组蛋白）

中文名称：ORM1 抗原（重组蛋白）

英文名称：ORM1 Antigen (Recombinant Protein)

别名：orosomuroid 1; ORM; AGP1; AGP-A

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 19-201 amino acids of human ORM1

技术规格：

Full name:	orosomuroid 1
Synonyms:	ORM; AGP1; AGP-A
Swissprot:	P02763
Gene Accession:	BC026238
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression.