

ORM2 抗原（重组蛋白）

中文名称：ORM2 抗原（重组蛋白）

英文名称：ORM2 Antigen (Recombinant Protein)

别名：AGP2; AGP-B; AGP-B'

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 19-201 amino acids of human ORM2

技术规格：

Full name:	orosomucoïd 2
Synonyms:	AGP2; AGP-B; AGP-B'
Swissprot:	P19652
Gene Accession:	BC015964
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression.