

PIP 抗原（重组蛋白）

中文名称： PIP 抗原（重组蛋白）

英文名称： PIP Antigen (Recombinant Protein)

别名： GPIP4; GCDFP15; GCDFP-15

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 29-146 amino acids of human PI

技术规格：

Full name:	prolactin-induced protein
Synonyms:	GPIP4; GCDFP15; GCDFP-15
Swissprot:	P12273
Gene Accession:	BC010950
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The hormonally responsive prolactin-inducible protein (PIP) gene is expressed in benign and malignant breast tumor tissues and in some normal exocrine organs such as sweat, salivary, and lacrimal glands.