

兔抗 EVPL 多克隆抗体

中文名称：兔抗 EVPL 多克隆抗体

英文名称：Anti-EVPL rabbit polyclonal antibody

抗 原：EVPL

别 名：envoplakin; EVPK

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the plakin family of proteins that forms a component of desmosomes and the epidermal cornified envelope. This gene is located in the tylosis oesophageal cancer locus on chromosome 17q25, and its deletion is associated with both familial and sporadic forms of oesophageal squamous cell carcinoma. Patients suffering from the autoimmune mucocutaneous disorder, paraneoplastic pemphigus, develop antibodies against the encoded protein.
Applications:	ELISA, IHC
Name of antibody:	EVPL

Immunogen:	Synthetic peptide of human EVPL
Full name:	envoplakin
Synonyms:	EVPK
SwissProt:	Q92817
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

