

兔抗 EVX1 多克隆抗体

中文名称: 兔抗 EVX1 多克隆抗体

英文名称: Anti-EVX1 rabbit polyclonal antibody

别名: even-skipped homeobox 1; EVX-1

抗原: EVX1

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the even-skipped homeobox family characterized by the presence of a homeo domain closely related to the Drosophila even-skipped (eve) segmentation gene of the pair-rule class. The encoded protein may play an important role as a transcriptional repressor during embryogenesis.
Applications:	ELISA, IHC
Name of antibody:	EVX1
Immunogen:	Synthetic peptide of human EVX1
Full name:	even-skipped homeobox 1
Synonyms:	EVX-1

SwissProt:	P49640
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human breast cancer
IHC Recommend dilution:	20-100



