

## RBBP8 抗原(重组蛋白)

中文名称: RBBP8 抗原(重组蛋白)

英文名称: RBBP8 Antigen (Recombinant Protein)

别名: RIM; COM1; CTIP; JWDS; SAE2; SCKL2

储 存: 冷冻(-20℃)

相关类别: 抗原

概 述:

Fusion protein corresponding to a region derived from 603-902 amino acids of human RBBP8

## 技术规格:

Full name:	retinoblastoma binding protein 8				
Synonyms:	RIM; COM1; CTIP; JWDS; SAE2; SCKL2				
Swissprot:	Q99708				
Gene Accession:	BC030590				
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE				
Expression system:	Escherichia coli				
Tags:	His tag C-Terminus, GST tag N-Terminus				
Background:	The protein encoded by this gene is a ubiquitously expressed nucl ear protein. It is found among several proteins that bind directly t o retinoblastoma protein, which regulates cell proliferation. This protein complexes with transcriptional co-repressor CTBP. It is also a ssociated with BRCA1 and is thought to modulate the functions of BRCA1 in transcriptional regulation, DNA repair, and/or cell cycle checkpoint control. It is suggested that this gene may itself be a t				



umor	suppressor	acting	in	tha	cama	nathway	20	RPC A1
ulliol	3uppressor	acting	1111	uie	Sairie	patriway	as	DIVCAT.