

兔抗 EXO5 多克隆抗体

中文名称：兔抗 EXO5 多克隆抗体

英文名称： Anti-EXO5 rabbit polyclonal antibody

别名： exonuclease 5; DEM1; Exo V; hExo5; C1orf176

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： EXO5

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a single-stranded DNA (ssDNA)-specific exonuclease that can slide along the DNA before cutting it. However, human replication protein A binds ssDNA and restricts sliding of the encoded protein, providing a 5'-directionality to the enzyme. This protein localizes to nuclear repair loci after DNA damage.
Applications:	ELISA, IHC
Name of antibody:	EXO5
Immunogen:	Fusion protein of human EXO5
Full name:	exonuclease 5
Synonyms:	DEM1; Exo V; hExo5; C1orf176

SwissProt:	Q9H790
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human cervical cancer
IHC Recommend dilution:	25-100



