

SCGB1A1 抗原（重组蛋白）

中文名称：SCGB1A1 抗原（重组蛋白）

英文名称： SCGB1A1 Antigen (Recombinant Protein)

别名： secretoglobin, family 1A, member 1 (uteroglobin); UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 22-91 amino acids of human SCGB1A1

技术规格：

Full name:	secretoglobin, family 1A, member 1 (uteroglobin)
Synonyms:	UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP
Swissprot:	P11684
Gene Accession:	BC004481
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.