

兔抗 EXOSC3 多克隆抗体

中文名称：兔抗 EXOSC3 多克隆抗体

英文名称： Anti-EXOSC3 rabbit polyclonal antibody

别名： exosome component 3; p10; PCH1B; RRP40; Rrp40p; CGI-102; hRrp-40; bA3J10.7

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： EXOSC3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a non-catalytic component of the human exosome, a complex with 3'-5' exoribonuclease activity that plays a role in numerous RNA processing and degradation activities. Related pseudogenes of this gene are found on chromosome 19 and 21. Alternatively spliced transcript variants encoding different isoforms have been described.
Applications:	ELISA, IHC
Name of antibody:	EXOSC3
Immunogen:	Fusion protein of human EXOSC3
Full name:	exosome component 3

Synonyms:	p10; PCH1B; RRP40; Rrp40p; CGI-102; hRrp-40; bA3J10.7
SwissProt:	Q9NQT5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human cervical cancer
IHC Recommend dilution:	25-100



