

## 兔抗 EXT1 多克隆抗体

- 中文名称：兔抗 EXT1 多克隆抗体
- 英文名称：Anti-EXT1 rabbit polyclonal antibody
- 别名：exostosin glycosyltransferase 1; EXT; LGS; TTV; LGCR; TRPS2
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：EXT1
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes an endoplasmic reticulum-resident type II transmembrane glycosyltransferase involved in the chain elongation step of heparan sulfate biosynthesis. Mutations in this gene cause the type I form of multiple exostoses.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	EXT1
<b>Immunogen:</b>	Fusion protein of human EXT1
<b>Full name:</b>	exostosin glycosyltransferase 1
<b>Synonyms:</b>	EXT; LGS; TTV; LGCR; TRPS2
<b>SwissProt:</b>	Q16394

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	50-200

