

兔抗 EXTL1 多克隆抗体

中文名称：兔抗 EXTL1 多克隆抗体

英文名称：Anti-EXTL1 rabbit polyclonal antibody

别名：EXTL

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：EXTL1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the multiple exostoses (EXT) family of glycosyltransferases, which function in the chain polymerization of heparan sulfate and heparin. The encoded protein harbors alpha 1,4- N-acetylglucosaminyltransferase activity, and is involved in chain elongation of heparan sulfate and possibly heparin.
Applications:	ELISA, IHC
Name of antibody:	EXTL1
Immunogen:	Synthetic peptide of human EXTL1
Full name:	exostosin-like glycosyltransferase 1
Synonyms :	EXTL

SwissProt:	Q92935
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

