

兔抗 EZH2 多克隆抗体

中文名称：兔抗 EZH2 多克隆抗体

英文名称： Anti-EZH2 rabbit polyclonal antibody

别名： WWS; ENX1; EZH1; KMT6; WWS2; ENX-1; EZH2b; KMT6A

相关类别： 一抗

抗原： EZH2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric protein complexes, which are involved in maintaining the transcriptional repressive state of genes over successive cell generations. This protein associates with the embryonic ectoderm development protein, the VAV1 oncoprotein, and the X-linked nuclear protein. This protein may play a role in the hematopoietic and central nervous systems. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene.

Applications:	ELISA, WB
Name of antibody:	EZH2
Immunogen:	Fusion protein of human EZH2
Full name:	enhancer of zeste 2 polycomb repressive complex 2 subunit
Synonyms:	WWS; ENX1; KMT6; WWS2; ENX-1; EZH2b; KMT6A
SwissProt:	Q15910
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	85 kDa
WB Positive control:	Raji and K562 cells
WB Recommended dilution:	200-1000

