

## SIGLEC9 抗原（重组蛋白）

中文名称： SIGLEC9 抗原（重组蛋白）

英文名称： SIGLEC9 Antigen (Recombinant Protein)

别名： CD329; CDw329; FOAP-9; siglec-9; OBBP-LIKE

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 146-336 amino acids of human SIGLEC9

技术规格：

<b>Full name:</b>	sialic acid binding Ig-like lectin 9
<b>Synonyms:</b>	CD329; CDw329; FOAP-9; siglec-9; OBBP-LIKE
<b>Swissprot:</b>	Q9Y336
<b>Gene Accession:</b>	BC035365
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Sialic acid-binding Ig-like lectin 9 is a protein that in humans is encoded by the SIGLEC9 gene. Putative adhesion molecule that mediates sialic-acid dependent binding to cells. Preferentially binds to alpha-2,3- or alpha-2,6-linked sialic acid. The sialic acid recognition site may be masked by cis interactions with sialic acids on the same cell surface. Expressed by peripheral blood leukocytes (neutrophils and monocytes).

ocytes but not eosinophils).