

兔抗 EZR(Ab-478) 多克隆抗体

中文名称：兔抗 EZR(Ab-478) 多克隆抗体

英文名称： Anti-EZR(Ab-478) rabbit polyclonal antibody

别名： CVL; CVIL; VIL2; HEL-S-105

抗原： EZR(Ab-478)

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The cytoplasmic peripheral membrane protein encoded by this gene functions as a protein-tyrosine kinase substrate in microvilli. As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this gene.
Applications:	WB

Name of antibody:	EZR(Ab-478)
Immunogen:	Synthesized non-phosphopeptide derived from human Ezrin around the phosphorylation site of tyrosine 478 (P-V-Y(p)-E-P).
Full name:	ezrin
Synonyms :	CVL; CVIL; VIL2; HEL-S-105
SwissProt:	P15311
WB Predicted band size:	69 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000

