

兔抗 EZR (Phospho-Thr567)多克隆抗体

- 中文名称: 兔抗 EZR (Phospho-Thr567)多克隆抗体
- 英文名称: Anti-EZR (Phospho-Thr567) rabbit polyclonal antibody
- 别名: CVL; CVIL; VIL2
- 相关类别: 一抗
- 储存: 冷冻 (-20℃) 避光
- 宿主: Rabbit
- 抗原: EZR (Phospho-Thr567)
- 反应种属: Human, Mouse, Rat
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	Probably involved in connections of major cytoskeletal structures to the plasma membrane. In epithelial cells, required for the formation of microvilli and membrane ruffles on the apical pole. Along with PLEKHG6, required for normal macropinocytosis.
Applications:	WB
Name of antibody:	EZR (Phospho-Thr567)
Immunogen:	Synthetic peptide of human EZR (Phospho-Thr567)
Full name:	ezrin (Phospho-Thr567)
Synonyms :	CVL; CVIL; VIL2
SwissProt:	P15311

WB Predicted band size:	69 kDa
WB Positive control:	MDA cells untreated or treated with anisomycin
WB Recommended dilution:	500-1000

