

## SPATA16 抗原（重组蛋白）

中文名称：SPATA16 抗原（重组蛋白）

英文名称： SPATA16 Antigen (Recombinant Protein)

别名： spermatogenesis associated 16; SPGF6; NYD-SP12

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 320-569 amino acids of human SPATA16

技术规格：

<b>Full name:</b>	spermatogenesis associated 16
<b>Synonyms:</b>	SPGF6; NYD-SP12
<b>Swissprot:</b>	Q9BXB7
<b>Gene Accession:</b>	BC034496
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a testis-specific protein that belongs to the tetratricopeptide repeat-like superfamily. The encoded protein localizes to the Golgi apparatus and may play a role in spermatogenesis.