

## 兔抗 F5(Coagulation factor V heavy chain) 多克隆抗体

中文名称：兔抗 F5(Coagulation factor V heavy chain)多克隆抗体

英文名称：Anti-F5(Coagulation factor V heavy chain) rabbit polyclonal antibody

别名：coagulation factor V; FVL; PCCF; THPH2; RPRGL1

抗原：F5(Coagulation factor V heavy chain)

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene encodes an essential cofactor of the blood coagulation cascade. This factor circulates in plasma, and is converted to the active form by the release of the activation peptide by thrombin during coagulation. This generates a heavy chain and a light chain which are held together by calcium ions. The activated protein is a cofactor that participates with activated coagulation factor X to activate prothrombin to thrombin. Defects in this gene result in either

	an autosomal recessive hemorrhagic diathesis or an autosomal dominant form of thrombophilia, which is known as activated protein C resistance.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	F5(Coagulation factor V heavy chain)
<b>Immunogen:</b>	Synthetic peptide of human F5(Coagulation factor V heavy chain)
<b>Full name:</b>	coagulation factor V
<b>Synonyms:</b>	FVL; PCCF; THPH2; RPRGL1
<b>SwissProt:</b>	P12259
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	20-100

