

## 兔抗 FAAH 多克隆抗体

中文名称：兔抗 FAAH 多克隆抗体

英文名称： Anti-FAAH rabbit polyclonal antibody

别名： PSAB; FAAH-1

相关类别： 一抗

抗原： FAAH

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that is responsible for the hydrolysis of a number of primary and secondary fatty acid amides, including the neuromodulatory compounds anandamide and oleamide.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	FAAH
<b>Immunogen:</b>	Fusion protein of human FAAH
<b>Full name:</b>	fatty acid amide hydrolase
<b>Synonyms :</b>	PSAB; FAAH-1
<b>SwissProt:</b>	O00519
<b>ELISA Recommended dilution:</b>	2000-5000

<b>IHC positive control:</b>	Human liver cancer and Human brain
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	63 kDa
<b>WB Positive control:</b>	LoVo cells, human testis tissue and A549 cells
<b>WB Recommended dilution:</b>	500-2000



