

兔抗 FAAP24 多克隆抗体

中文名称：兔抗 FAAP24 多克隆抗体

英文名称：Anti-FAAP24 rabbit polyclonal antibody

别名：Fanconi anemia core complex associated protein 24; C19orf40

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FAAP24

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Plays a role in DNA repair through recruitment of the FA core complex to damaged DNA. Regulates FANCD2 monoubiquitination upon DNA damage. Induces chromosomal instability as well as hypersensitivity to DNA cross-linking agents, when repressed. Targets FANCM/FAAP24 complex to the DNA, preferentially to single strand DNA.
Applications:	ELISA, IHC
Name of antibody:	FAAP24
Immunogen:	Synthetic peptide of human FAAP24
Full name:	Fanconi anemia core complex associated protein

	24
Synonyms:	C19orf40
SwissProt:	Q9BTP7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human lung cancer
IHC Recommend dilution:	30-150



