

TAX1BP1 抗原（重组蛋白）

中文名称：TAX1BP1 抗原（重组蛋白）

英文名称：TAX1BP1 Antigen (Recombinant Protein)

别名：Tax1 (human T-cell leukemia virus type I) binding protein 1; T6BP; TXBP151; CALCOCO3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 180-420 amino acids of human TAX1BP1

技术规格：

Full name:	Tax1 (human T-cell leukemia virus type I) binding protein 1
Synonyms:	T6BP; TXBP151; CALCOCO3
Swissprot:	Q86VP1
Gene Accession:	BC050358
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a HTLV-1 tax1 binding protein. The encoded protein interacts with TNFAIP3, and inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. Degradation of this protein by caspase-3-like family proteins is associated with apoptosis induced by TNF. This protein may also have a role in the inhibition of inflammatory signaling pathways. Alternatively spliced transcript varia

nts encoding different isoforms have been found for this gene.[]