

TEKT1 抗原（重组蛋白）

中文名称：TEKT1 抗原（重组蛋白）

英文名称：TEKT1 Antigen (Recombinant Protein)

别名：Tektin 1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 300 amino acids of human TEKT1

技术规格：

Full name:	Tektin 1
Swissprot:	Q969V4
Gene Accession:	BC014599
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene product belongs to the tektin family of proteins. Tektns comprise a family of filament-forming proteins that are coassembled with tubulins to form ciliary and flagellar microtubules. This gene is predominantly expressed in the testis and in mouse, tektin 1 mRNA was localized to the spermatocytes and round spermatids in the seminiferous tubules, indicating that it may play a role in spermatogenesis.