

## TLR8 抗原（重组蛋白）

中文名称：TLR8 抗原（重组蛋白）

英文名称：TLR8 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概 述：

Fusion protein corresponding to C terminal 200 amino acids of human TLR8

技术规格：

<b>Full name:</b>	toll-like receptor 8
<b>Synonyms:</b>	CD288
<b>Swissprot:</b>	Q9NR97
<b>Gene Accession:</b>	BC101075
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a member of the Toll-like receptor (TLR) family which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns (PAMPs) that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. The various TLRs exhibit different patterns of expression. This gene is predominantly expressed in lung and peripheral blood leukocytes, and lies in close proximity

to another family member, TLR7, on chromosome X.