

TNFAIP1 抗原（重组蛋白）

中文名称：TNFAIP1 抗原（重组蛋白）

英文名称：TNFAIP1 Antigen (Recombinant Protein)

别名：tumor necrosis factor, alpha-induced protein 1 (endothelial); B12; B61; EDP1; BTBD34; hBACURD2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 32-316 amino acids of human TNFAIP1

技术规格：

Full name:	tumor necrosis factor, alpha-induced protein 1 (endothelial)
Synonyms:	B12; B61; EDP1; BTBD34; hBACURD2
Swissprot:	Q13829
Gene Accession:	BC001643
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene was identified as a gene whose expression can be induced by the tumor necrosis factor alpha (TNF) in umbilical vein endothelial cells. Studies of a similar gene in mouse suggest that the expression of this gene is developmentally regulated in a tissue-specific manner.

mlbio 梅联生物
Good elisakit producers