

兔抗 FAIM 多克隆抗体

中文名称：兔抗 FAIM 多克隆抗体

英文名称：Anti-FAIM rabbit polyclonal antibody

别名：FAIM1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FAIM

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene protects against death receptor-triggered apoptosis and regulates B-cell signaling and differentiation. Several transcript variants encoding different isoforms have been found for this gene
Applications:	ELISA, WB, IHC
Name of antibody:	FAIM
Immunogen:	Synthetic peptide of human FAIM
Full name:	Fas apoptotic inhibitory molecule
Synonyms :	FAIM1
SwissProt:	Q9NVQ4

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human colon cancer and Human cervical cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	20 kDa
WB Positive control:	Raji cells
WB Recommended dilution:	200-1000





