

TNFSF13 抗原（重组蛋白）

中文名称：TNFSF13 抗原（重组蛋白）

英文名称：TNFSF13 Antigen (Recombinant Protein)

别名：APRIL; CD256; TALL2; ZTNF2; TALL-2; TRDL-1; UNQ383/PRO715

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 105-247 amino acids of human TNFSF13

技术规格：

Full name:	tumor necrosis factor (ligand) superfamily, member 13
Synonyms:	APRIL; CD256; TALL2; ZTNF2; TALL-2; TRDL-1; UNQ383/PRO715
Swissprot:	O75888
Gene Accession:	BC008042
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the tumor necrosis factor (TNF) ligand family. This protein is a ligand for TNFRSF17/BCMA, a member of the TNF receptor family. This protein and its receptor are both found to be important for B cell development. In vitro experiments suggest

ted that this protein may be able to induce apoptosis through its interaction with other TNF receptor family proteins such as TNFRSF6/FAS and TNFRSF14/HVEM. Alternative splicing results in multiple transcript variants. Some transcripts that skip the last exon of the upstream gene (TNFSF12) and continue into the second exon of this gene have been identified; such read-through transcripts are contained in GeneID 407977, TNFSF12-TNFSF13.