

## 兔抗 FAM20C 多克隆抗体

中文名称：兔抗 FAM20C 多克隆抗体

英文名称： Anti-FAM20C rabbit polyclonal antibody

别名： FAM20C, golgi associated secretory pathway kinase; RNS; DMP4; G-CK; DMP-4; GEF-CK

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： FAM20C

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the family of secreted protein kinases. The encoded protein binds calcium and phosphorylates proteins involved in bone mineralization. Mutations in this gene are associated with the autosomal recessive disorder Raine syndrome.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FAM20C
<b>Immunogen:</b>	Synthetic peptide of human FAM20C
<b>Full name:</b>	FAM20C, golgi associated secretory pathway kinase

<b>Synonyms:</b>	RNS; DMP4; G-CK; DMP-4; GEF-CK
<b>SwissProt:</b>	Q8IXL6
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-200

