

TRAF2 抗原（重组蛋白）

中文名称：TRAF2 抗原（重组蛋白）

英文名称：TRAF2 Antigen (Recombinant Protein)

别名：TNF receptor-associated factor 2; TRAP; TRAP3; MGC:45012

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 351-496 amino acids of human TRAF2

技术规格：

Full name:	TNF receptor-associated factor 2
Synonyms:	TRAP; TRAP3; MGC:45012
Swissprot:	Q12933
Gene Accession:	BC032410
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the TNF receptor associated factor (TRAF) protein family. TRAF proteins associate with, and mediate the signal transduction from members of the TNF receptor superfamily. This protein directly interacts with TNF receptors, and forms a heterodimeric complex with TRAF1. This protein is required for TNF-alpha-mediated activation of MAPK8/JNK and NF-kappaB. The protein complex formed by this protein and TRAF1 interacts with the inhibitor-of-apoptosis proteins (IAPs), and functions as a me

diator of the anti-apoptotic signals from TNF receptors.