

## 兔抗 FAM35A 多克隆抗体

中文名称: 兔抗 FAM35A 多克隆抗体

英文名称: Anti-FAM35A rabbit polyclonal antibody

别名: FAM35A1; bA163M19.1

抗原: FAM35A

储存: 冷冻 (-20°C) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	Belongs to the FAM35 family. 2 isoforms of the human protein are produced by alternative splicing.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	FAM35A
<b>Immunogen:</b>	Synthesized peptide derived from internal of human FAM35A/B.
<b>Full name:</b>	family with sequence similarity 35, member A
<b>Synonyms :</b>	FAM35A1; bA163M19.1
<b>SwissProt:</b>	Q86V20/A6NDK8
<b>WB Predicted band size:</b>	94 kDa
<b>WB Positive control:</b>	Jurkat cells lysate

**WB Recommended dilution:**

500-3000

