

兔抗 FAM50A 多克隆抗体

中文名称：兔抗 FAM50A 多克隆抗体

英文名称：Anti-FAM50A rabbit polyclonal antibody

别名：family with sequence similarity 50 member A; 9F; XAP5; HXC26; HXC-26; DXS9928E

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FAM50A

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the FAM50 family. The encoded protein is highly conserved in length and sequence across different species. It is a basic protein containing a nuclear localization signal, and may function as a DNA-binding protein or a transcriptional factor.
Applications:	ELISA, WB, IHC
Name of antibody:	FAM50A
Immunogen:	Synthetic peptide of human FAM50A
Full name:	family with sequence similarity 50 member A
Synonyms:	9F; XAP5; HXC26; HXC-26; DXS9928E

SwissProt:	Q14320
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	40 kDa
WB Positive control:	Hela cell lysate
WB Recommended dilution:	500-2000



