

TRAF3IP3 抗原（重组蛋白）

中文名称：TRAF3IP3 抗原（重组蛋白）

英文名称：TRAF3IP3 Antigen (Recombinant Protein)

别名：TRAF3 interacting protein 3; T3JAM

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 240-436 amino acids of human TRAF3IP3

技术规格：

Full name:	TRAF3 interacting protein 3
Synonyms:	T3JAM
Swissprot:	Q9Y228
Gene Accession:	BC130417
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The gene encodes a protein that mediates cell growth by modulating the c-Jun N-terminal kinase signal transduction pathway. The encoded protein may also interact with a large multi-protein assembly containing the phosphatase 2A catalytic subunit. Alternative splicing results in multiple transcript variants.

mlbio 梅联生物
Good elisakit producers